Smilyss®

Ready to Smile

by **SPORT-ELEC**
INSTITUT

FR

ES

NL

DE

GB

RO

Mode d'emploi

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Smilyss®

SPORT ELEC® Institut - 27520 BOURGTHEROULDE - FRANCE

www.sport-elec.com

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IPX7

IMPORTANT SAFETY INSTRUCTIONS

Before using your toothbrush, please pay particular attention to the instruction manual, it contains important information on the product use, safety, care and maintenance.

- Please store this instruction manual in a safe place so you can refer to the manual if you need to.
- It is advisable to supply the device to operated circuit-breakers without integral and a rated power releasing no more than 30mA.
- Only use this device with alternating current 230V. Before plugin, please check the voltage of your home matches the voltage device. You can refer to the device identification sheet (Charger: AC 230 volts, alternating current (~), 50/60HZ., 2.5 W – Power Pack: Ni-MH 1/2AA 300mAh 2.4V). The device must be plugged in a safe standards electrical distribution connection.

CAUTIONS:

Never immerse the appliance in water or any other liquid.

Do not use abrasive cleaners or household solvents.

- If the appliance stops unexpectedly during the use, failed or if the device or its cord are damaged, switch off and unplug the device from the outlet.
 - Always unplug the device from the electrical outlet before cleaning or maintenance, and after use.
- Do not let appliance cords hang from counter-tops.
- Keep the appliance and cable away from sources of heat.
- Always use a suited and all working extension cord.
- Do not pull the cord to remove the appliance from the outlet.
- Keep the appliance away and out of reach of children.
- Do not use in wet environments (sauna, hammam) and near electrical appliance or trunk main with earthing connection.
- Never immerse the appliance in water or any other liquid.
- Do not use outdoors. This appliance has been designed for domestic indoor use only.
- Do not put the device down when it is in motion.
- Only use the accessories provided with the appliance.
- Switch off the power supply and unplug from the wall before cleaning and or changing accessories.
- Do not operate your appliance in the presence of explosive or smoke.
- Any intervention or repair must be done by an authorized repairer: unauthorized repairer may represent a real danger for the user.

GB

- This appliance has been designed for domestic indoor use only. Inadequate or improper use of the device will cancel the warranty and won't commit the liability of the manufacturer.
- The supply cord cannot be replaced. If the cord is damaged the appliance should be scrapped.

CONNECTING

The cord of the base can be plug in any alternating current 220-240 volts (~) including the shaver outlets. Place the toothbrush on the charging base. When charging the red light indicator goes on. The power pack may become hot when charging. It is not affecting its functioning, overload is excluded. In order to keep your toothbrush power full; always put it back on its base.

MAINTENANCE

Before using the appliance for the first time and to protect the life of the battery, we recommend to proceed as mentioned below:

1. Charge the battery entirely for at least 16 hours prior to use.
2. We recommend that you run the battery flat by putting on the toothbrush until it stops automatically. The toothbrush has a timer for 2 minutes. To run the battery flat entirely you need to launch the program several times.
3. Then fully charge the battery for at least 16 hours.

This procedure "charging/recharging" must be done each 6 months.

The duration of use on full power is around 20-25 minutes.

N.B. when the battery starts running flat, the toothbrush stops automatically to avoid batteries to be fully discharged. All you need is plug the device back on the base to recharge.

INSTRUCTIONS MANUAL

- After operating the "charging /recharging" step mentioned in the previous paragraph.
- Screw the brush head on the power pack.
- Before the first use of a brush head, carefully rinse the brush head.
- Rinse your mouth before brushing.
- To avoid spillage: Place the brush at an angle to your teeth, then put the toothbrush on (button on the handle of the toothbrush while you're at a lower speed-the first blue LED light).

ENVIRONMENT FRIENDLY DISPOSAL

Your appliance contains valuable materials which can be recovered or recycled.
Please leave it at a collection point as listed by your local authority.

CE (2004/108/CE)-(2006/95/CE)

Rohs (2002/95/EC)

IPX7

- Pass the brush on your teeth from top to bottom. Guide it along the inner and outer faces also taking care to clean the gums. It is not necessary to press hard on the brush. Clean the sides and surfaces of your teeth with a twisting motion.
- If you want to upgrade to a higher speed pressing of the button on the handle of the toothbrush (turbo speed - 2 blue LED are lit) and then to go to massage gums press once again (a single blue LED lit top of the brush).
- To stop the toothbrush again press the button (all LED off).
- To fully remove dental plaque, the timer allows you to brush your teeth for 2 min.

N.B. *clinical tests shown that the efficient way to clean your teeth is to gently run the brush over from one tooth to another, without pressing.*

- Carefully rinse your mouth after brushing.
- During the first use, your gums may bleed a little particularly if they are sensitive. After few days the bleeding should stop. If your gums keep bleeding after 15 days, we recommend you to consult your dentist.
The shape of the brush is changing with use, which depends on the personal way of use of the brush. This misshapeness has no impact on the brush efficiency. However do not forget to change the brush head as soon as a manual toothbrush.
- After using the toothbrush, rinse carefully the head brush.
- Put the toothbrush back on the charging base.

CARE AND CLEANING

● HEAD BRUSH:

To keep your brush head in a good condition of use, rinse it carefully under water during few seconds.

The toothbrush has to be in motion while operating.

For a perfect hygiene, you can soak the brush head (and only the head) once a week in a water container with a denture cleaner.

● TOOTHBRUSH:

Clean the toothbrush with a soft and damp cloth.

The power pack has to be perfectly dry before putting the toothbrush back on the charging base.

● CHARGING BASE:

Clean the external sides from the charging base with a soft and damp cloth.

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